

STRIA DRAFT PROCEDURES

Test ID: BLANK
STRIA

Saturday 1st January, 2022

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden


SUP: David

OSO: David/Remko

TC: Kira

TO: Cristina

CP: Lukas



Authors:

Lukas Welzel, Kira Strelow, Cristina Cordun

Version:

V1.0

In case of emergency: +31 (0)71 527 5737

Low risk	Medium risk	High risk
<ul style="list-style-type: none"> • Test-setup is safe to approach • No safety gear required 	<ul style="list-style-type: none"> • Only authorized personnel in test area • Wear appropriate safety gear 	<ul style="list-style-type: none"> • Clear all personnel from test area • Do not approach the test-setup



Abbreviations

CAM	Camera / livestream person	HE	Helium gas
CP	Command post	LN2	Liquid Nitrogen
OSO	Operations Safety Officer	LOX	Liquid oxygen
SUP	External Supervisor	MA	Mirror Assembly
TC	Test Conductor	N2	Nitrogen gas
TL	Test Leader	KE	Knife Edge
TO	Test Operator	OP	Optical Platform
		RTFM	Read The Full Manual
		TS	Translation Stage

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA
Test ID: BLANK

Summary of the setup

The physical system is a schlieren imager consisting of two separate assemblies: the optical platform and the mirror assembly.

By cutting the incoming beam in half, differences in the refractive index over some area of the media cross section in the beam are shown as brighter or darker regions.

In these procedures the assembly, calibration and usage of the STRIA system are covered.

Test Summary

The main goal of this test is to image a medium in flux:

1. Assembly:
2. Preparation: Pack the STRIA system and the sample.
3. Calibration: Calibrate the STRIA system for the specific test setup.
4. Test: Image the sample.
5. Clean-up.

Test Goals and Operations

The following primary and secondary goals are as follows:

Primary Goals:

- Image flow of a medium
- Find [TBD for specific test]
- Verify safe operations of the STRIA system
- Verify the timing sequence

Secondary Goals:

- Verify STRIA procedures
- Have fun

The following tests are preliminary covered by these procedures:

- TBD: Compressed N2 Up, Slow
- TBD: Compressed N2 Up, Fast
- TBD: Compressed N2 Down, Slow
- TBD: Compressed N2 Down, Fast
- TBD: Candle
- TBD: Flow over Surface
- TBD: Soldering Iron
- TBD: Flow over Cylinder

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

Automated Sequence

The automated sequence that will be run has the following timeline:

Time [ms]	Action
-500	Start high speed data logging
0	Release the medium
Test Duration	Stop releasing the medium
Test Duration + 500	Stop high speed data logging
Attention	Timing preliminary, to be determined with test data

Emergency and Unpowered System States

When an emergency button is pressed the systems will go to the emergency state. When all power is cut, they will go to their unpowered state.

System	Emergency state	Unpowered state
CAM	NA	Off
LED	NA	Off
Arduino	NA	Off
Laptop	NA	Off

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

System Schematics

Below a diagram of the OP is shown.




Figure 1: Schematic of the OP.

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

“We can easily forgive a child who is afraid of the dark; the real tragedy of life is when men are afraid of the light.”

Plato

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

Packing list

A - General				
ID	Check		Description	Comments
A1			PR kit	
A2			Nitrile gloves	
A3			at least 3 x Clipboards	
A4			at least 3 x Procedures	print single sided, two copies (TC, SO) need to be printed in color

B - Test Setup				
ID	Check		Description	Comments
B1			MA	
B2			OP	

C - Electronics				
ID	Check		Description	Comments
C1			Laptop + charger	
C2			USB - Ethernet adapter	
C3			HDD/USB stick/SD card with min 64 GB	only needed if the data needs the be recovered from the laptop immediately

D - Safety				
ID	Check		Description	Comments
D1			First Aid Kit	
D2			≥ 2 sets of nitrile gloves	
D3			CO2 Fire extinguisher	

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA
Test ID: BLANK

Checklist

E - Preparation, day before leaving			
ID	Check		Description
E1			Confirm all items on packing list are packed
E2			Check again

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

Assembly

Item Lists for Assembly

F - Item List for MA				
ID	Check		Description	Comments
F1			Mirror	
F2			2 x Mirror Mount Leg (Thorlabs)	
F3			Mirror Mount	
F4			4 x M6x20	mounting the MA to the optical table
F5			4 x M6x16	mounting the mirror mount to the mirror mount legs
F6			2 x M6x70 (partially threaded)	mounting the mirror mount to the mirror mount legs

Lego pieces are given as WxLxH.

G - Item List for OP				
ID	Check		Description	Comments
G1			200x200 Lego platform	
G2			25 x 2x4x3 Lego pieces	
G3			25 x 2x4x1 Lego pieces	
G4			25 x 1x4x3 Lego pieces	
G5			25 x 1x4x1 Lego pieces	
G6			25 x 2x2x3 Lego pieces	
G7			25 x 2x2x1 Lego pieces	
G8			6 x 6x30x1 Lego pieces	
G9			2x single axis LEGO TS	
G10			1x two axis LEGO TS	
G11			2x KE base	
G12			2x KE holder	
G13			4x M3x30	
G14			2x M3 spring	
G15			4x M3 washer	
Continued				

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden


Saturday 1st January, 2022

STRIA

Test ID: BLANK

ID	Check	Description	Comments
G16		2x Lens mount base	
G17		1x Lens mount	
G18		1x camera mount	
G19		1x lens	TODO: right lens name
G20		4x rubber bands	
G21		1x KE	
G22		1x LED	
G23		4x M6x80	
G24		4x M6 medium spring	
G25		4x M6 washer	
G26		4x M6 heat set threaded inserts	
G27		5x M3 heat set threaded inserts	
G28		LED assembly	

Assembly Procedures

H - MA Assembly			
ID	Check	Description	Comments
		Do not touch optical surfaces. Only touch optics using gloves but do not touch their optical surfaces. RTFM of the optics.	
			
H1		Bolt the long side of the legs to the mirror mount	
H2		Bolt the short side of the legs to the optical table	
H3		Insert the mirror carefully into the mirror mount	you will need at least 2 people for the following steps
H4		Make sure the mirror sits flat against the mirror mount	
H5		Fasten the mirror using the plastic screws	
H6		Fasten the protective screen to the mirror mount	

In case of emergency: +31 (0)71 527 5737



Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

I - OP Assembly				
ID	Check	Description	Comments	
<div></div>				
		The heat set threaded inserts will be pressed in using a soldering iron. RTFM of the inserts.		
<div></div>				
I1			Press the M6 heat set threaded inserts into the two lens bases from the top using a soldering iron	make sure the top of the inserts sit flush against the countersunk holes of the base
I2			Press the M3 heat set threaded inserts into the two KE holder bases from the top using a soldering iron	make sure the top of the inserts sit flush against the countersunk holes of the base
I3			Press the M3 heat set threaded inserts into the camera holder from the top using a soldering iron	make sure the top of the inserts sit flush against the countersunk holes of the base
I4			Let the parts with the inserts cool down	
I5			Check the alignment of the inserts using bolts 1. Thread a long bolt into the insert 2. Check from all sides if the bolt is sufficiently straight (normal to the surface) 3. If the bolt is not sufficiently straight: (a) Remove the bolt (b) Heat up the insert and improve the alignment by pushing on it with the soldering iron (c) Let the insert cool down and check straightness (d) If aligning the insert is not possible in that way heat up the insert, quickly insert the bolt and straighten the bolt and insert	make sure the top of the inserts always sit flush against the base
I6			Push the M3 bolts through the KE and LED holder	
I7			Insert a M3 washer, a M3 spring and again a M3 washer onto the M3 bolts	
I8			thread the bolts into the KE and LED bases	
Continued				

In case of emergency: +31 (0)71 527 5737



Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

ID	Check	Description	Comments
I9		Push the M6 bolts through the bolt holes of the lens and camera holders	
I10		Insert a M6 washer, a medium M6 spring and again a M6 washer onto the M6 bolts	
I11		Thread the M6 bolts into the lens and camera base threaded inserts	
		Do not touch optical surfaces. Only touch optics using gloves but do not touch their optical surfaces. RTFM of the optics.	
 			
I12		Insert the lens into the lens mount	the strongly convex side should face away from the lens mount
I13		Insert the lens holder back-plate into the lens mount	be careful to not put too much pressure on the lens
I14		Secure the back-plate to the lens mount using rubber bands	Use at least 4 bands holding it in different places
I15		Carefully make sure the lens sits securely in the mount	
I16		Insert the camera into the camera mount	the c-mount should be aligned with the flat side of the mount
I17		Fasten the camera using three M3 bolts and the heat threaded inserts on top	you might need an insert to protect the camera from the bolts
I18		Mount the sub-assemblies on the Lego TS: <ul style="list-style-type: none">• KE: 1 axis (x-direction)• Lens mount: 1 axis (y-direction)• Camera mount: 2 axis (x, y direction)	
I19		Mount the sub-assemblies on the Lego base according to Figure 1	
I20		Mount the LED (-assembly) on the LED holder	
I21		Make sure everyone is aware of the KE being mounted on the OP by informing the TC, SO, SUP that the KE will now be mounted	
I22		Mount the KE on the KE holder	

In case of emergency: +31 (0)71 527 5737


Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

J - Calibration			
ID	Check	Description	Comments
			
System Status: <ul style="list-style-type: none"> • System is off • Camera is off • LED is off • Mirror cover is mounted to mount • Knife edge might be uncovered 			
J1		Notify CP that calibration will start	
J2		Be aware of the potentially exposed, sharp knife edge	
J3		Make sure that the area between the OP and MA is clear of obstacles	
J4		Make sure that the MA is securely mounted	
J5		Make sure that the cover of the camera is screwed into the C mount (detector is covered)	
J6		Make sure that the OP is securely mounted	
J7		Make sure that the OP is securely assembled	
J8		Make sure that the cable running to the camera is securely plugged into the camera and laptop	
J9		Make sure that the cables running from the LED to the breadboard are securely plugged in	
J10		Make sure that the cables running from the breadboard to the Arduino are securely plugged in	
J11		Make sure that the cable running from the Arduino to the laptop is plugged in	
J12		Position the OP at $2400 \text{ mm} \pm 5 \text{ mm}$ from the MA, measured from the OP-knife edge/LED to the MA front plane	
J13		Align the OP and MA by eye <ol style="list-style-type: none"> 1. The optical axis of the MA should go through the midpoint of the KE-LED line 2. The LED-MA-KE triangle should be symmetric about the MA optical axis 	
J14		Start the laptop	
Continued			

In case of emergency: +31 (0)71 527 5737


Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

ID	Check	Description	Comments
J15		Start the SpinView program	
J16		Start the camera	
J17		Connect the camera to SpinView	You should see a flat, black image
J18		Turn on the camera live view (do not start recording)	
J19		Remove the camera cover from the camera	SpinView should show a gray image
The LED will now be powered on			
			
J20		Power on the LED	
J21		Align the LED 1. Place a screen behind the MA 2. Turn the light in the lab off 3. Bend the LED until it is roughly centered on the center of the mirror 4. Turn the light in the laboratory back on	
J22		Adjust the distance of the OP to the MA until the spot of the reflected LED is on or near the KE	use a screen between the KE and lens if necessary
J23		Adjust the position of the KE (x, z direction) until the spot is focused on the KE and half of the light is blocked	
J24		Place a screen between the lens and the camera and check if a spot is visible	
The detector of the camera will now be exposed			
J25		Move the mirror on the z and y axis until the light beam is centered on it	
J26		Move the camera (y,z direction) until the image (mirror and/or sample) is shown	
J27		Place an object with sharp features in the optical path where the sample will be placed	
J28		Move the camera on the x direction until the image of the object becomes sharp	
J29		Remove the object from the optical path	
J30		If the image is still not sharp return to J26	
J31		Make a test recording using the python or Spin-View scripts	
Continued			

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA
Test ID: BLANK

ID	Check	Description	Comments
J32		Make sure that the recording is successful	
J33		If the image and video is sharp the system is calibrated	
J34		Turn the camera to standby (turn off live feed)	
J35		Turn off the LED	
J36		(Re)Mount the mirror cover to the MA	
System Status: <ul style="list-style-type: none"> • System is on standby • System is calibrated • Camera is on standby • LED is on • Mirror cover is mounted to MA 			

In case of emergency: +31 (0)71 527 5737


Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

K - Perform Imaging			
ID	Check	Description	Comments
			
System Status: <ul style="list-style-type: none"> • System is on standby • System is calibrated • Camera is on standby • LED is on • Mirror cover is mounted to MA 			
K1		Announce that imaging will be performed	
K2		Prepare the python script to record images	
K3		(optional) Remove the mirror cover	
K4		Ascertain that everybody stands still during the test and only the TO is close to the setup	
K5		Announce that testing will start now	
K6		Count down	
K7		Start high speed recording	
K8		Start flow/other test	
K9		Announce that the test is concluded	
K10		(optional) re-mount the mirror cover	
K11		Ascertain that the data has been saved	
System Status: <ul style="list-style-type: none"> • System is on standby • System is calibrated • Camera is on standby • LED is on • Mirror cover is mounted to MA 			

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
Niels Bohrweg 2
2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK

L - Clean up and Setup Disassembly.			
ID	Check	Description	Comments
Tentative procedures.			
L1		Verify that the system is safe to work on	
L2		Cover the knife edge and take it out of the setup	
L3		Cover the optical instruments	camera, lens
L4		Disconnect all electronics	
L5		Cover the setup	
L6		Clean up test area 1. No trash 2. No exposed optics 3. No loose bolts 4. No loose Lego	
System Status: <ul style="list-style-type: none"> • System is packed up 			

In case of emergency: +31 (0)71 527 5737

Test Location:

Optical Laboratory, Faculty of Science, Leiden University
 Niels Bohrweg 2
 2333CA Leiden

Saturday 1st January, 2022

STRIA

Test ID: BLANK